

The high security connection



# Preparation devices for medications/solutions

# Bag

Brand name	Code	Product	Description	Indication	Quantity box/case
qimo <b>spike</b> B	7218.01		qimo <b>spike</b> B is an IV fluid bag access spike fitted with qimo Q.	Aspiration or addition of medication/ fluid from/into an IV fluid bag.	50/600
qimo <b>spike</b> B valve	7218.11	Maria	qimo <b>spike</b> B valve is an IV fluid bag access spike fitted with qimo $\varphi$ and a non-return valve.	Aspiration only from an IV fluid bag, thus avoiding the accidental contamination of the IV solution.	25/450

#### Vial

Brand name	Code	Product	Description	Indication	Quantity box/case
qimo <b>cap</b> ø14	7219.01 7219.010*		qimo <b>cap</b> is a vial access cap fitted with qimo Q. It has been specifi-	Withdrawal and reconstitution of small volumes of	50/900 4/72
qimo <b>cap</b> ø20	7219.02 7219.020*		cally designed to create a permanent connection with the medication vial.	medication in a completely closed environment.	25/450 4/48
qimo <b>spike</b> V	7217.04 7217.040*		qimo <b>spike</b> V is a vented spike for vials fitted with qimo $\varphi$ .	Repeated with- drawals from every type of vial.	25/300 4/48
qimo <b>harpoon</b>	7213.02		qimoharpoon is a closed system transfer device that has an air-intake tube connected	Withdrawal and reconstitution	10/80
qimo <b>harpoon</b> S	7213.01		to a 50ml balloon filled with sterile air which ensures pressure equalization.	of high-risk medications.	10/00

## **Syringe/Extension tube**

Brand name	Code	Product	Description	Indication	Quantity box/case
qimo <b>3</b>	7210.91 7210.910*		qimo of is a disinfectable closed needleless safety connector. It is a safety captive adaptor for male Luer-lock hubs. Once attached, it cannot be disconnected. It connects only to qimo of the connects only to qimo of the closed is a disconnected in the connects only to qimo of the closed is a disconnected in the connects only to qimo of the closed is a disconnected in the closed in the clos	To secure all the male Luer-lock fittings.	50/900 4/72
qimo <b>transfer</b>	7220.01	***************************************	Female Luer-lock to female Luer-lock connector, with qimo Q on both ends of the connector.	To transfer medica- tions from syringe to syringe.	50/600

<sup>\*</sup>Multipack by 25

#### Elastomeric pump - 2 packs available

	Brand name	Code	Product	Description	Indication	Quantity box/case
Option 1	qimo <b>prime,</b> qimo & & qimo P	7210.04		Set with gimo <b>prime</b> , gimo $\delta$ and gimo $\Phi$ in the same packaging.	Can be used to secure the filling port of the elastomeric pump and the distal end of its tubing, but other applications are possible.	50/600
Option 2	qimo 🗣	7210.02		qimo Q is a multi- purpose, closed, needleless and disinfectable safety connector, supplied in individual packaging.	qimo Q can be used to connect the tubing of the distal end of the elastomeric pump secured with qimo o <sup>7</sup> , to access an implantable port or any other patient vascular device.	50/1200
	qimo <b>prime</b> & qimo <b>3</b>	7210.01		Set supplied with qimo of and qimo-prime in the same packaging.	Can be used to secure distal end of elastomeric pump tubing.  Need to be connected to a qimo Q	50/600





# Transfer devices for medications/solutions

### Bag

Brand name	Code	Product	Description	Indication	Quantity box/case
qimo <b>connect</b>	7221.01 7221.020*		IV fluid bag extension set with a lateral injection site fitted with qimo Q. qimo $\delta$ secures the distal end of the extension tube.	Dilution and administration of all kind of medications.	25/150 10/60
qimo <b>connect</b> filter	7221.011 7221.0110*		IV fluid bag extension set with a 0.2µm filter and a lateral injection site fitted with qimo Q. qimo $\delta$ secures the distal end of the extension tube.	Dilution and administration of solutions which need to be filtered (ex: Paclitaxel).	25/150 10/60

<sup>\*</sup>Multipack by 10 (without qimo 3)



### **Generic products**



Brand name	Code	Product	Description	Indication	Quantity box/case
qimo <b>pump</b>	7227.101 7227.201 7227.401		qimo <b>pump</b> is a safe adaptor for pump infusion sets. This device is fitted with 1, 2 or 4 qimo Q connectors, a spike with air-venting port (with a cap) and a universal tube to connect to pump dedicated sets.	To secure the dedicated pump infusion sets and connect the bag containing medication in a safe manner.	25/150
qimo <b>IVset</b> with integral qimo <b>octopus</b>	7229.101 7229.201 7229.401		Intravenous infusion device with spike and drip chamber. It is fitted with 1, 2 or 4 qimo Q and qimo do attach to the male Luerdal and the discount of the dis	cations by gravity or	25/150
qimo <b>IVset</b> without qimo <b>octopus</b>	7229.102 7229.202 7229.402	1011	lock port at the distal end of gimo <b>IVset</b> . gimo <b>prime</b> allows the priming of the line once gimo $\delta$ has been attached.	with qimono infusion pump (QIP).	
qimo <b>connect</b> 오	7221.02		IV fluid bag extension tube with a lateral injection site fitted with qimo Q and an air-venting port (with a cap).	nistration of non toxic	25/150
qimo <b>prime</b> & qimo <b>3</b>	7210.01		Set supplied with gimo of and gimoprime in the same packaging.		50/600
qimo <b>octopus</b>	7231.221		qimooctopus is a safe double-lumen IV extension tube allowing continuous infusion, intermittent injection and blood sampling. Its lateral line is fitted with qimo Q. qimo dhas to be attached to the distal end of qimooctopus. qimoprime allows the priming of the line once qimo dhas been attached.	Safe emergency access. Easy direct injection with a syringe close to the patient.	50/300







Brand name	Code	Product	Description	Indication	Quantity box/case
		Infusion set for g	gravity or peristaltic pump		
qimo <b>IVset</b> with integral qimo <b>octopus</b>	7229.101 7229.201 7229.401		Intravenous infusion device with spike and drip chamber. It is fitted with 1, 2 or 4 qimo Q and qimo 3 to attach to the	To infuse safely medications by gravity or	25/150
qimo <b>IVset</b> without qimooctopus	7229.102 7229.202 7229.402	10111	male Luer-lock port at the distal end of qimo <b>IVset</b> . qimo <b>prime</b> allows the priming of the line once qimo $\delta$ has been attached.	with gimono infusion pump (QIP).	23/130
		A	Accessory		
qimono infusion pump ( <b>QIP</b> )	0VEPQIM FR 0VEPQIMG02 EN PLUG 0VEPQIMC03 NL 0VEPQIMC04 DE 0VEPQIMC06 ES 0VEPQIMC07 PT 0VEPQIMC08 IT 0VEPQIMC10 PL	© more than the state of the st	Infusion Pump: qimono infusion pump (QIP) is safe, simple to use, precise, silent. 2 set up modes available: Volume & Time of infusion or Flow rate & Volume to infuse. It displays a large coloured screen, visual + audible critical alarms for nurse monitoring & patient safety. It is easy to clean and offers a quick access to the Lithium battery. 2 years warranty (excluding battery).	Qimono Infusion Pump is dedicated to administration of cytotoxic drugs and perfusion in day care hospital.  To be used ONLY in conjonction with qimo IVset listed above	1 unit/box

#### Patient's vascular access



Brand name			Code		Product	Description	Indication	Quantity box/case
qimo <b>perfusafe</b>			7242.720 7242.725 7242.920 7242.925			qimoperfusafe is a safety infusion set with Huber needle fitted with a qimo Q to access implantable ports.	Safety device (protected needle and closed system connector) to puncture an implantable port.	10/90
qimo P			7210.02			qimo Q is a multi-purpose, closed, needleless and disinfectable safety connector.	To connect to any other vascular access, Huber needles, CVCsafely and with a closed system	50/1200
	Gauge	Ø	Longueur	Code				
			15mm	741507*				
			17mm	741707				
	22G	0,7	20mm	742007		Safety huber needle with		
			25mm	742507				
			30mm	743007			qimo <b>flow</b> +	
			35mm	743507		extension tube and safety	safety Huber	
			15mm	741509*		adaptor: qimo	needles are indicated	
			17mm 20mm	741709 742009		우 (closed system)	for the	40.65
qimo <b>flow</b> +*	20 <b>G</b>	0,9	25mm	742009		qimo <b>flow</b> + is	administration or withdrawal	12/30
			30mm	743009		a single-handed activation	of fluids	
			35mm	743509		needle allowing	through implantable	
			15mm	741511*		an automatic positive flush	ports.	
			17mm	741711		during removal		
	466	4.4	20mm	742011				
	19G	1,1	25mm	742511				
			30mm	743011				
			35mm	743511				

 $<sup>\</sup>ensuremath{^{*}}$  Please check availability with your Vygon local representative





#### **Vesical instillation**

Brand name	Code	Product	Description	Indication	Quantity box/case
qimo <b>adapt</b>	7210.06		qimo <b>adapt</b> is a safe adaptor for the connection of a syringe on the funnel of a tube.  qimo Q allows to connect a syringe secured with qimo o <sup>7</sup> .	Safe bladder instillation	50/600

# **Priming device**

Brand name	Code	Product	Description	Indication	Quantity box/case
qimo <b>prime</b>	7210.81		qimo <b>prime</b> is a priming device which can be connected to any product fitted with qimo $\sigma^{-1}$ .	Priming device for qimo o <sup>71</sup>	50/1200

#### **Premedication extension line**

Brand name	Code	Product	Description	Indication	Quantity box/case
tubing set	855.02		Tubing set is an extension tube with vented spike and a male Luer-lock connector with a non-return valve and a priming cap.	Extension tube for drug administration	25/150











Notes		
	_	
	_	
	_	
	_	
	_	
	-	
	_	







#### For further information, please contact: questions@vygon.com

The specifications shown in this leaflet are for information only and are not, under any circumstances, of a contractual nature.

Vygon – 5, rue Adeline • 95440 ECOUEN • FRANCE Reception: +33 (0)1.39.92.63.63 - Service clients France: +33 (0)1.39.92.63.81 Export customer service: +33 (0)1.39.92.64.15 Fax.: +33 (0)1.39.92.64.44 • www.vygon.com









