## SWISS DUCCLAST® ADIAL & FOCUSED, THE ULTIMATE SHOCK WAVE THERAPY SOLUTION





# SWISS DUOCLAST®

> Always the right technology at hand for tendinopathies, muscle pain or bone pathologies: shallow with Radial Shock Waves and deep with Focused Shock Waves.



#### **SWISS DOLORCLAST® MASTER** → ALL SUPERFICIAL INDICATIONS

**SWISS PIEZOCLAST®** → ALL DEEP INDICATIONS

 $\mathbf{O}$ . 0



#### BALLISTIC GENERATION

BADIAL ESW1

> Compressed air accelerates a projectile which strikes a fixed applicator. The kinetic energy is converted into a Shock Wave delivered to the target tissue through the skin.

#### PIEZOCERAMIC GENERATION

> High voltage is applied to 1,000 piezoceramic crystals generating pressure waves, which penetrate the skin and travel through the tissue focusing on a cigar-shaped volume.

### SWISSDOLORCLAST.COM



For more information: welcome@ems-ch.com